


Glucose Urine Timed

Order Name: **GLUC TM U**
Test Number: 2053700
Revision Date: 05/20/2021

TEST NAME	METHODOLOGY	LOINC CODE
Total Urine Volume		3167-4
Creatinine, Urine, mg/dL	Kinetic Alkaline Picrate (Jaffe)	20624-3
Creatinine Urine Timed Estimated		
Creatinine, Urine, 24 Hour		2162-6
Glucose 24 Hour Urine mg/dL	Hexokinase	21305-8
Estimated Glucose 24 Urine		2351-5
Glucose 24 Hour Urine mg/hr		2350-7
Hours Collected - Prompt		30211-7
Volume 24 Urine Prompt		

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Urine, 24-hour	24 hour Urine Container	Refrigerated
Instructions	Timed urine collection. Record number of hours and volume in ml on the specimen container. Urine specimens should be maintained at refrigerated temperature and analyzed as soon as possible. Stability: Ambient 2 hours, Refrigerated 2 hours, Frozen 2 days			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Used with non diabetic patients whose urine tests positive for reducing substances
CPT Code(s)	82945 (81050 if urine measured by RML)
Service Provided By	 labcorp Oklahoma, Inc.