

Chlorpromazine (Thorazine) Level

Order Name: **THORAZINE**
 Test Number: 2008650
 Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Chlorpromazine (Thorazine) Level	<u>Liquid Chromatography/Tandem Mass Spectrometry</u>	3471-0

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (0.5 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Alternate 1	2 mL (0.5 mL)	Plasma	Lithium Heparin (Dark Green Top / No-Gel)	Refrigerated
Alternate 2	2 mL (0.5 mL)	Plasma	EDTA (Lavender Top)	Refrigerated
Instructions	<p>Specimen Type: Red-top tube, lavender-top (EDTA) tube OR green-top (heparin) tube. DO NOT USE A GEL-BARRIER TUBE. The use of gel-barrier tubes is not recommended due to slow absorption of the drug by the gel. Depending on the specimen volume and storage time, the decrease in drug level due to absorption may be clinically significant.</p> <p>Specimen Storage: Room Temperature</p> <p>Specimen Collection: Transfer separated serum or plasma to a plastic transport tube.</p> <p>Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION

Expected TAT	3 - 5 days
Notes	Labcorp Test Code: 072132
CPT Code(s)	80342

Service Provided By  Oklahoma, Inc.