

**Cortisol Free, 24-Hour Urine**

Order Name: **CORT FR U**  
Test Number: 3602275  
Revision Date: 03/01/2023

TEST NAME	METHODOLOGY	LOINC CODE
Cortisol Free, 24-Hour Urine	<a href="#">Liquid Chromatography/Mass Spectroscopy</a>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	100 mL aliquot (10 mL aliquot )	Urine, 24-hour	24 hour Urine Container	Room Temperature
Instructions	<p><b>Notes:</b> 10 mL aliquot (Note: This volume Does NOT allow for repeat testing.)</p> <p><b>Specimen Type:</b> Plastic urine container</p> <p><b>Container Detail:</b> No Preserv</p> <p><b>Specimen Storage:</b> Room Temperature</p> <p><b>Specimen Collection:</b> Instruct the patient to void at 8 AM and discard the specimen. Then collect all urine including the final specimen voided at the end of the 24-hour collection period (ie, 8 AM the next morning) in a plastic bottle. Screw the lid on securely. (It is acceptable to collect cortisol with boric acid or HCl preservative.) HCl preservative or without a preservative.)</p> <p><b>Special Instructions:</b> State 24-hour volume.</p> <p><b>Specimen Stability:</b> Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION	
Expected TAT	4 - 6 days
Notes	Labcorp Test Code: 004432
CPT Code(s)	82530
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.