St. John Health System
Lab Catalog

Catecholamines, Fractionated, Plasma

Order Name: CAT FRAC P
Test Number: 3801400
Revision Date: 12/10/2022

TEST NAME		ME'	THODOLOGY	LOINC CODE
Catecholamines, Fractionated, Plasma		Hig	h Performance Liquid Chromatography	49257-9
SPECIMEN REQUIREM	MENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (2.2 mL)	Plasma	EDTA (Lavender Top)	Frozen
Instructions	Specimen Type: Lavender-top (EDTA) tube or green-top (heparin) tube Specimen Storage: FREEZE. After centrifugation, the plasma can be stored up to two hours at room temperature. Sample can be kept for up to two weeks at -20(degrees)C. Specimen Collection: Draw blood in lavender-top (EDTA) tube or green-top (heparin) tube. Invert to mix with preservatives. Centrifuge and transfer the plasma to labeled plastic transport tube. Freeze immediately (within one hour after collection) at -20(degrees)C and ship frozen. The time between blood collection and the preparation of plasma is quite critical; if the time exceeds one hour, catecholamine values increase (when blood is kept at 4(degrees)C) or decrease (when left at 20(degrees)C).[1] To avoid delays in turnaround time when requesting multiple tests on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH TEST REQUESTED.			

GENERAL INFORMATION		
Expected TAT	4-6 Days	
Notes	Labcorp Test Code: 084152	
CPT Code(s)	82384	
Service Provided By	labcorp Oklahoma, Inc.	

Specimen Stability: Ambient: Not Available, Refrigerated: Not Available, Frozen: Not Available