

Serotonin 5-HIAA, 24-Hour (5-Hydroxyindoleacetic Acid)


Order Name: **SEROTON U**
Test Number: 3808450
Revision Date: 03/17/2023

TEST NAME	METHODOLOGY	LOINC CODE
Serotonin 5-HIAA, 24-Hour (5-Hydroxyindoleacetic Acid)	High Performance Liquid Chromatography	

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	4 mL (1 mL)	Urine, 24-hour	24 hour Urine Container	Room Temperature
Instructions	<p>This is not for Random Urine Collections</p> <p>Notes: 1 mL aliquot (Note: This volume Does NOT allow for repeat testing).</p> <p>Specimen Type: Plastic urine container, no preservative (NOTE: 1 g/L boric acid may be added as a preservative for other tests without harm to 5-HIAA.)</p> <p>Container Detail: No preserv</p> <p>Specimen Storage: Room Temperature</p> <p>Specimen Collection: Instruct the patient to void at 8 AM and discard the specimen. Then collect all urine including the final specimen voided at the end of the 24-hour collection period (ie, 8 AM the next morning). Screw the lid on securely. Label container with patient's name, date, and time. Measure and record total urine volume. Mix well; send aliquot.</p> <p>Special Instructions: Measure and record total 24-hour urine volume on the request form.</p> <p>Specimen Stability: Ambient: 7 days, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION

Expected TAT	4 - 8 days
Notes	Labcorp Test Code: 004069
CPT Code(s)	83497
Service Provided By	 labcorp Oklahoma, Inc.