

Hepatitis C Genotype

Order Name: **HCV GENO**Test Number: 5594650
Revision Date: 01/14/2025

| TEST NAME | | METHODOLOGY | | LOINC CODE |
|----------------------|--|---------------------------|---------------------|-----------------------|
| Hepatitis C Genotype | | Polymerase Chain Reaction | | 32286-7 |
| SPECIMEN REQUIRE | MENTS | | | |
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 2.5 mL (1.5 mL) | Plasma | EDTA (Lavender Top) | Frozen |
| Alternate 1 | 2.5 mL (1.5 mL) | Serum | Clot Activator SST | Frozen |
| Instructions | Preferred Specimen is 2.5 mL (1.5 mL) Plasma or Serum from EDTA (Lavender Top) or Clot Activator SST Centrifuge specimen and separate plasma or serum from cells within 6 hours of collection. Serum or plasma specimens may be stored frozen at -20°C or lower for up to 60 days. Avoid multiple freeze/thaws. | | | |

| GENERAL INFORMATION | | |
|---------------------|---|--|
| Testing Schedule | Thursday | |
| Expected TAT | may not be reported until next day | |
| Clinical Use | HCV viral genotype is used to predict the likelihood of therapeutic response and determine duration of treatment. Patient must have a viral load of > 500 IU/mL for Genotype to be determined. | |
| CPT Code(s) | 87902 | |
| Service Provided By | Oklahoma, Inc. | |