

Occult Blood X1

Order Name: **OCC BL 1**
Test Number: 3510000
Revision Date: 01/10/2019

TEST NAME	METHODOLOGY	LOINC CODE
Occult Blood X1	<u>Guaiac Colormetric Reaction (GUIAC)</u>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.5 mL (0.1)	Stool, Random	Fecal Occult Blood Card	Room Temperature
Instructions	<p>SPECIAL DIET REQUIRED:</p> <p>Do not eat red meat, any blood-containing food, cantaloupe, uncooked broccoli, turnip, radish, or horseradish for 3 days prior to the test. Some medications may interfere with this test. These include vitamin C and aspirin. The health care provider should be consulted regarding medication changes that may be necessary. Medication should not be stopped or decreased without consulting the health care provider.</p> <p>The patient's full name and date/time of collection should be noted on the card.</p> <p>Stability: Specimen should be immediately transferred to the provided sample card by the patient or provider at time of specimen collection. Occult blood cards should be submitted to the lab within 10 days of the first sample collection. To insure stability, submit cards promptly after the last collection. Cards containing sample can be stored at room temperature for up to 14 days before testing.</p>			

GENERAL INFORMATION	
Testing Schedule	Mon-Fri
Expected TAT	1-3 Days
Clinical Use	Colon cancer screening
CPT Code(s)	82270
Service Provided By	 labcorp Oklahoma, Inc.