

Hepatitis C RNA Quantitative PCR with reflex to Genotyping

Order Name: **HCV Q+GENO**Test Number: 5594677
Revision Date: 09/10/2019

| TEST NAME | METHODOLOGY | LOINC CODE |
|--------------------------------|---------------------------|------------|
| HCV PCR Quantitative IU/mL | Polymerase Chain Reaction | 11011-4 |
| HCV PCR Quantitative Log IU/mL | Polymerase Chain Reaction | 38180-6 |
| | | |

| SPECIMEN REQUIREMENTS | | | | | |
|---|--|------------------------------|---|-----------------------|--|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment | |
| Preferred | 2.5 mL (1.5 mL) | Plasma | EDTA (Lavender Top) | Frozen | |
| Alternate 1 | 2.5 mL (1.5 mL) | Serum | Clot Activator SST | Frozen | |
| Instructions | For HCV Genotype the viral load must be at least 500 IU/mL | | | | |
| | Preferred Specimen is 2.5 m | L (1.5 mL) Plasma or Serum f | rom EDTA (Lavender Top) or Clot Activator | SST | |
| Centrifuge specimen and separate plasma or serum from cells within 2 hours of collection. | | | | | |
| | Serum or plasma specimens may be stored frozen at -20°C or lower for up to 60 days. Avoid multiple freeze/thaws. | | | /thaws. | |
| | | | | | |

| GENERAL INFORMATION | |
|---------------------|--|
| Testing Schedule | Thursday |
| Expected TAT | may not be reported until next day |
| Clinical Use | HCV RNA viral genotype is used to predict the likelihood of therapeutic response and determine duration of treatment. HCV genotype testing will be performed only if an HCV viral load is detected via PCR testing. The linear range is 12 to 100,000,000 IU/mL (1.08 to 8.00 Log IU/mL). Order Code: [5594650] Hepatitis C Genotyping Result Code: [5594650] HCV Geno |
| Notes | If the HCV genotype cannot be determined by PCR, the specimen will be sent out for HCV genotype sequencing analysis. An additional CPT code 87902 will apply. |
| CPT Code(s) | 87522 (w/ possible 87902 for genotype) |
| Service Provided By | Oklahoma, Inc. |