

Heavy Metals Panel, Urine (24hr or Random) with Reflex to Arsenic Fractions

Specimen Stability: Ambient: 14 days, Refrigerated: 14 days, Frozen: 14 days

Order Name: METAL S U
Test Number: 3810575
Revision Date: 11/15/2023

TEST NAME			METHODOLOGY	LOINC CODE	
Heavy Metals Panel, Urine (24hr or Random) with Reflex to Arsenic Fractions			Atomic Absorption Spectrophotometry		
SPECIMEN REQUIRE	EMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	15 mL (5.3 mL)	Urine, 24-hour	24 Hour Urine Acid Washed Container	Room Temperature	
Alternate 1	15 mL (5.3 mL)	Urine, Random	Acid Washed, Trace Element Free Contatiner	Room Temperature	
Instructions	Notes: 5.3 mL (Note: This volume Does NOT allow for repeat testing.)				
	Specimen Type: Plastic urine container, no preservative				
	Specimen Storage: Maintain specimen at room temperature.				
	Specimen Collection: Sampling time is not critical for industrial monitoring. OPTIONAL PROTOCOL: Instruct the patient to void at 8 AM and discard				
	the specimen. Then collect all urine including the final specimen voided at the end of the 24-hour collection period (ie, 8 AM the next morning). Avoid				
	contact with metal during collection. Screw the lid on securely. See individual tests for collection times for industrial monitoring.				
	Special Instructions: Reques	t form MUST state 24-hour s	pecimen volume, if applicable.		

St. John Health System Lab Catalog

GENERAL INFORMATION

Expected TAT

3 - 6 days

Notes

If total arsenic concentration is between 35-2000 ug/L, then Arsenic, Fractionated shown below, will be added to determine the proportion of organic, inorganic, and methylated forms. Additional charges apply.

Arsenic, Fractionated [Order Code: 3709000]

RESULT CODE	RESULT NAME	LOINC
3800475	Arsenic Inorganic	12481-8
3800276	Arsenic, MMA	94458-7
3800476	Arsenic, DMA	94457-9
3800676	Arsenic, Total Toxic	n/a
3800876	Arsenic, Please Note:	n/a
3800875	Arsenic Organic	n/a
3800675	Arsenic Methylated	n/a

CPT Code(s)

82570, 83825, 82175, 83655

Service Provided By

