St. John Health System Lab Catalog

Fragile X Syndrome, Carrier

Order Name: **FRAGILE X**Test Number: 0117075
Revision Date: 03/11/2024

TEST NAME		METHO	DDOLOGY	LOINC CODE	
Fragile X Syndrome, Carrier		PCR/C	apillary Electrophoresis	53039-4	
SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	8.5 mL (3 mL)	Whole Blood	EDTA (Lavender Top)	Room Temperature	
Alternate 1	1	See Instructions	See Instructions	Room Temperature	
Instructions	Specimen Type: 8.5 mL whole blood from Lavender-top (EDTA) tube OR PurFlock buccal swab kit or Oragene Dx saliva kit  Collection: Standard phlebotomy. Follow PurFlock buccal swab kit or Oragene Dx 500 saliva kit collection instructions. Do not eat, drink, smoke or chew gum 30 minutes prior to collection.  Specimen Storage: Room temperature  Stability: Whole blood: 14 days at room temperature or 30 days at 4'C				

GENERAL INFORMATION		
Expected TAT	6-14 Days	
Clinical Use	This test is used for carrier screening for fragile X syndrome (FMR1).  ?Limitations: Technologies used do not detect germline mosaicism and do not rule out the presence of large chromosomal aberrations including rearrangements and gene fusions, or variants in regions or genes not included in this test, or possible inter/intragenic interactions between variants, or repeat expansions. Variant classification and/or interpretation may change with time if more information becomes available. False positive or false negative results may occur for reasons that include: rare genetic variants, sex chromosome abnormalities, pseudogene interference, blood transfusions, bone marrow transplantation, somatic or tissue-specific mosaicism, mislabeled samples or erroneous representation of family relationships. This test was developed and its performance characteristics determined by Labcorp. It has not been cleared or approved by the Food and Drug Administration.	
Notes	Polymerase chain reaction (PCR) followed by capillary electrophoresis, with reflex to AGG interruption analysis and methylation PCR analysis for positive samples	
CPT Code(s)	81243	
Service Provided By	Oklahoma, Inc.	