St. John Health System Lab Catalog

Hepatitis C Analyzer Reflex to Genotype

Order Name: **HEP C AN+G**Test Number: 5590857
Revision Date: 05/28/2021

TEST NAME		N	IETHODOLOGY	LOINC CODE
Hepatitis C Analyzer Reflex to Genotype			Chemiluminescence Assay	
SPECIMEN REQUIR	EMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	See Instructions	Serum and Plasma	See Instructions	See Instructions
Instructions	Two Specimens are Required for this Analyzer.			
	Specimen 1 - for the Hepatitis C Antibody 4mL(2mL) Serum from Clot Activator SST, Keep Refrigerated or Plasma from Lithium Heparin PST (Light Green) Stability: Ambient 3 days, refrigerated 7 days, Frozen >7 days. Specimen 2 - for possible HCV Quantitative Viral Load by PCR			(Light Green)
	Centrifuge specimen and sep Alternately specimens may be	parate plasma or serum from ce	m EDTA (Lavender Top) or Clot Active ells within 2 hours of collection. rs prior to separation and centrifugation. ours at 2-8°C.	

Serum or plasma specimens may be stored frozen at -20°C or lower for up to 60 days. Avoid multiple freeze/thaws.

GENERAL INFORMATION		
Testing Schedule	Mon-Fri	
Expected TAT	1-2 Days (Assay Dependant)	
Clinical Use	Positive Hepatitis C Antibody will order the following test for viral load and genotype determination. Hepatitis C RNA Quantitative PCR with reflex to Genotyping [5594677] HCV Q+GENO	
Notes	Assay Schedule: HCV Quantitative: Dayshift- Set up on Monday, Wednesday and Friday reported same day. HCV Genotype: Dayshift- Set up on Thursday and reported on Friday.	
CPT Code(s)	86803 (possible 87522 for HCV Quantitative Viral Load by PCR) (possible 87902 for HCV Genotyping)	
Service Provided By	labcorp Oklahoma, Inc.	