St. John Health System Lab Catalog

Gamma-Hydroxybutyric Acid (GHB)

Order Name: **GHB Serum**Test Number: 3604705
Revision Date: 10/01/2023

TEST NAME			METHODOLOGY	LOINC CODE
Gamma-Hydroxybutyric Acid (GHB)		Liquid Chromatography/Tandem Mass Spectrome	etry 59681-7	
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	5 mL (1.2 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Alternate 1	5 mL (1.2 mL)	Plasma	Lithium Heparin (Dark Green Top / No-Gel)	Refrigerated
Instructions	Specimen Type: Red-top tube or green-top (heparin) tube Specimen Storage: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Specimen Collection: Serum or plasma should be separated from cells within two hours of venipuncture. Send serum or plasma in a plastic transport tube. Gel-barrier tubes are NOT recommended. Special Instructions: Testing referred to MedTox Laboratories Inc TC 19250			

GENERAL INFORMATION	
Expected TAT	7 - 10 days
CPT Code(s)	80307
Service Provided By	labcorp Oklahoma, Inc.