

Tocainide Level

Instructions

Order Name: **TOCAINIDE**Test Number: 3611150
Revision Date: 12/12/2022

TEST NAME		METH	ODOLOGY	LOINC CODE
Tocainide Level		Liquid Chromatography/Tandem Mass Spectrometry		
SPECIMEN REQUIRE	MENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5 mL)	Serum	Clot Activator (Red Top, No-Gel)	Room Temperature
Alternate 1	1 mL (0.5 mL)	Plasma	Lithium Heparin (Dark Green Top / No-Gel)	Room Temperature
Alternate 2	1 mL (0.5 mL)	Plasma	EDTA (Lavender Top)	Room Temperature

Specimen Type: Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended.

Specimen Storage: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.

Specimen Collection: Serum or plasma should be separated from cells within two hours of venipuncture. Submit serum or plasma in a plastic transport tube.

Special Instructions: Testing referred to MedTox Laboratories Inc tc 121

GENERAL INFORMATION	
Expected TAT	5 - 7 days
Notes	Labcorp Test Code: 811372
CPT Code(s)	80299
Service Provided By	labcorp Oklahoma, Inc.