


**Propafenone Level**

Order Name: **PROPAFENON**  
Test Number: 3657050  
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Propafenone Level	<u>Liquid Chromatography/Tandem Mass Spectrometry</u>	6905-4

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1.0mL ( 0.5mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Alternate 1	1.0mL ( 0.5mL)	Plasma	Sodium Heparin (Green Top, No-Gel)	Refrigerated
Instructions	<p><b>Specimen Type:</b> Red-top tube or green-top (heparin) tube. Gel-barrier tubes are NOT recommended.</p> <p><b>Specimen Collection:</b> <b>Serum or plasma should be separated from cells within two hours of venipuncture.</b> Send serum or plasma in a plastic transport tube.</p> <p><b>Special Instructions:</b> Testing referred to MedTox Laboratories Inc TC 353</p> <p><b>Specimen Storage / Stability:</b> SUBMISSION/TRANSPORT (&lt;3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Specimen is stable for three days at room temperature, for seven days refrigerated, and for six months frozen.</p>			

GENERAL INFORMATION	
Expected TAT	7 - 14 days
Notes	Labcorp Test Code: 270088
CPT Code(s)	80299
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.