

**Acetazolamide Level**

Order Name: **Acetazolamide**  
Test Number: 4503200  
Revision Date: 12/20/2022

TEST NAME	METHODOLOGY	LOINC CODE
Acetazolamide Level	<u>Quantitative High Performance Liquid Chromatography</u>	3300-1

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1mL (0.5 mL)	Serum	Clot Activator (Red Top, No-Gel)	Room Temperature
Instructions	<p><b>Specimen Type:</b> Red-top tube or green-top (heparin) tube. Gel-barrier tubes are NOT recommended.</p> <p><b>Specimen Storage:</b> SUBMISSION/TRANSPORT (&lt;3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.</p> <p><b>Specimen Collection:</b> Serum or plasma should be separated from cells within two hours of venipuncture. Send serum or plasma in a plastic transport tube.</p> <p><b>Special Instructions:</b> Testing referred to MedTox Laboratories Inc MX#359</p>			

GENERAL INFORMATION	
Expected TAT	7-14 days
Notes	Labcorp Test Code: 811984
CPT Code(s)	80299
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.