


**Fluphenazine, Serum or Plasma**

Order Name: **Fluphenazine Lvl**  
Test Number: 6905015  
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Fluphenazine, Serum or Plasma	<u>Quant Liquid Chromatography-Tandem Mass Spectrometry</u>	3650-9

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (1 mL)	Serum	Clot Activator (Red Top, No-Gel)	Room Temperature
Alternate 1	2 mL (1 mL)	Plasma	Lithium Heparin PST (Light Green Top)	Room Temperature
Instructions	<p><b>Timing of specimen collection:</b> Pre-dose (trough) draw - At steady state concentration.</p> <p><b>Specimen Type:</b> Red-top tube, lavender-top (EDTA) tube OR green-top (heparin) tube. DO NOT USE A GEL-BARRIER TUBE. The use of gel-barrier tubes is not recommended due to slow absorption of the drug by the gel. Depending on the specimen volume and storage time, the decrease in drug level due to absorption may be clinically significant.</p> <p><b>Specimen Storage:</b> Room Temperature</p> <p><b>Specimen Collection:</b> Transfer serum or plasma to a transport tube with plastic stopper within 45 minutes of collection.</p> <p><b>Specimen Stability:</b> Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION	
Expected TAT	3 - 4 days
Notes	Labcorp Test Code: 706887
CPT Code(s)	80342
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.