

**Posaconazole Quantl by LC-MS/MS**

Order Name: **POSACON AF**  
Test Number: 6905121  
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Posaconazole Quantl by LC-MS/MS	Quantitative HPLC/Tandem Mass Spectrometry	53731-6

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.6 mL)	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Alternate 1	1 mL (0.6 mL)	Plasma	EDTA (Lavender Top)	Frozen
Instructions	<p><b>Specimen Type:</b> Red-top tube OR lavender-top (EDTA) tube</p> <p><b>Specimen Collection:</b> Transfer separated serum or plasma to a plastic transport tube. <b>DO NOT USE A GEL BARRIER TUBE.</b> The use of gel-barrier tubes is not recommended due to slow absorption of the drug by the gel. Depending on the specimen volume and storage times, the decrease in drug level due to absorption may be clinically significant.</p> <p><b>Specimen Storage / Stability:</b> Refrigerate Stability: Room temperature up to 14 day Refrigerate up to 14 days Frozen up to 14 days</p>			

GENERAL INFORMATION	
Expected TAT	2 - 7 days the next day.
Notes	Labcorp Test Code: 700265
CPT Code(s)	80187
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.