

Pyruvic Acid, CSF

Order Name: PYRU CSF
Test Number: 6905131
Revision Date: 12/12/2022

TEST NAME		METH	IODOLOGY	LOINC CODE	
Pyruvic Acid, CSF		Quantitative Enzymatic 14122-6		14122-6	
SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	1 mL (0.5)	CSF (Cerebrospinal Fluid)	Sterile Screwtop Container	Frozen	
Instructions	Specimen Type: Sterile screw capped vial  Specimen Storage: Freeze as soon as possibe after collection  Specimen Collection: Collect entire sample into a single sterile tube. Immediately after collection, add a 2:1 ratio (acid:sample) of COLD 8% perchloric acid. Vortex for 30 seconds, and place on ice for 5 minutes. Store frozen at -20 degrees Celsius and ship frozen.  Special Instructions: Testing referred to Medical Neurogenetic LLC MNEGA#MET11  Specimen Stability: Ambient: unstable (stability provided by manufacturer or li, Refrigerated : 24 hours (stability provided by manufacturer or li, Frozen: -20 C = 72 hours; -80 C = Indefinate (stability pr				

GENERAL INFORMATION		
Expected TAT	10 - 14 days	
Notes	Labcorp Test Code: 620045	
CPT Code(s)	84210	
Service Provided By	labcorp Oklahoma, Inc.	