

Tiagabine (Gabitril), Serum/Plasma

Order Name: **Tiagabine**
Test Number: **6905167**
Revision Date: **12/20/2022**

TEST NAME	METHODOLOGY	LOINC CODE
Tiagabine (Gabitril), Serum/Plasma	<u>Quantitative Liquid Chromatography/Tandem Mass Spectrometry</u>	21565-7

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	3 mL (0.6 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	<p>Specimen Type: Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended.</p> <p>Specimen Storage: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Specimen is stable for two weeks at room temperature, for two weeks refrigerated, and for seven months frozen.</p> <p>Specimen Collection: Serum or plasma should be separated from cells within two hours of venipuncture. Submit serum or plasma in a plastic transport tube.</p> <p>Special Instructions: Testing referred to MedTox Laboratories Inc MX#409</p> <p>Specimen Stability: Ambient: 2 weeks, Refrigerated : 2 weeks, Frozen: 7 months</p>			

GENERAL INFORMATION	
Expected TAT	4-7 days
Notes	Labcorp Test Code: 842104
CPT Code(s)	80199
Service Provided By	 labcorp Oklahoma, Inc.