

**Sertraline, Serum or Plasma**

Order Name: **Sertraline**  
Test Number: 7001007  
Revision Date: 11/14/2022

TEST NAME	METHODOLOGY	LOINC CODE
Sertraline, Serum or Plasma	<u>Quantitative Liquid Chromatography/Tandem Mass Spectrometry</u>	6906-2

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (0.5 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Instructions	<p><b>Specimen Type:</b> Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended.</p> <p><b>Specimen Storage:</b> SUBMISSION/TRANSPORT (&lt;3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.</p> <p><b>Specimen Collection:</b> Serum or plasma should be separated from cells within two hours of venipuncture. Submit serum or plasma in a plastic transport tube.</p> <p><b>Special Instructions:</b> Testing referred to MedTox Laboratories Inc MX#122</p>			

GENERAL INFORMATION	
Expected TAT	5-10 days
Notes	Labcorp Test Code: 808494
CPT Code(s)	80299
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.