

**Clobazam Quantitative, Serum or Plasma**

Order Name: **Clobazam Quant S/P**  
Test Number: 6906653  
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Clobazam	Quant Equilibrium Dialysis/HPLC-Tandem Mass Spectrometry	
N-Desmethyloclobazam	Quant Equilibrium Dialysis/HPLC-Tandem Mass Spectrometry	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Alternate 1	1 mL (0.5 mL)	Plasma	Lithium Heparin (Dark Green Top / No-Gel)	Refrigerated
Instructions	<p><b>Specimen Type:</b> Red-top tube or green-top (heparin) tube. <b>Gel-barrier tubes are NOT recommended.</b></p> <p><b>Specimen Storage:</b> SUBMISSION/TRANSPORT (&lt;3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Specimen is stable for nine days at room temperature, refrigerated, and frozen.</p> <p><b>Specimen Collection:</b> Serum or plasma should be separated from cells within two hours of venipuncture. Send serum or plasma in a plastic transport tube.</p> <p><b>Special Instructions:</b> Testing referred to MedTox Laboratories Inc TC 2837</p>			

GENERAL INFORMATION	
Expected TAT	4 - 10 days
Notes	Labcorp Test Code: 790500
CPT Code(s)	80339 (Alt code: G0480)
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.