

Eosinophil Cationic Protein

Order Name: **Eos Cation Prot**
Test Number: 6906929
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Eosinophil Cationic Protein	ImmunoCAP *****	25638-8

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.25 mL)	Serum	Clot Activator SST	Refrigerated or Frozen
Instructions	<p>Notes: 0.3 mL (Note: This volume Does NOT allow for repeat testing)</p> <p>Specimen Type: NMR lipoprofile serum-separator tube (LabCorp item ID 60360). Other tubes should not be used as they may produce erroneous results.</p> <p>Specimen Storage: Refrigerate</p> <p>Specimen Collection: Allow the sample to stand for 60 to 120 minutes to release the ECP in the cells prior to centrifugation. Separate the serum after 60 to 120 minutes and transfer to a plastic transport tube.</p> <p>Specimen Stability: Ambient: 5 days, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION	
Expected TAT	4 - 7 days
Notes	Labcorp Test Code: 004180
CPT Code(s)	83520
Service Provided By	 labcorp Oklahoma, Inc.