


**Lipoprotein Phenotyping**

Order Name: **LP PHENO**  
Test Number: 5194353  
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Lipoprotein Phenotyping	<a href="#">See Test Notes</a>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (1 mL)	Serum	Clot Activator (SST or Red No-Gel)	Refrigerated
Instructions	<p><b>Specimen Type:</b> Red-top tube or gel-barrier tube</p> <p><b>Specimen Storage:</b> Refrigerate</p> <p><b>Specimen Collection:</b> Separate serum from cells within 45 minutes of collection. Lipid panels are best avoided following acute myocardial infarct, for up to three months, although cholesterol can be measured in the first 24 hours.</p> <p><b>Special Instructions:</b> State patient's age, sex, and pertinent clinical history on the test request form.</p> <p><b>Specimen Stability:</b> Ambient: Unstable, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION	
Expected TAT	5 - 9 days
Notes	Electrophoretic with Frederickson phenotyping criteria Labcorp Test Code: 235036
CPT Code(s)	80061, 83700
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.