St. John Health System Lab Catalog

Bacterial Vaginosis, NAA

Order Name: Bact Vag NAA

Test Number: 5195322 Revision Date: 09/12/2023

TEST NAME METHODOLOGY LOINC CODE

Bacterial Vaginosis, NAA

Nucleic acid amplification (NAA)

SPECIMEN REQUIREMENTS Specimen Specimen Volume (min) Specimen Type Specimen Container Transport Environment Preferred 1 Swab APTIMA Multitest Swab -or- Unisex Swab Kit Room Temperature

Instructions Container: Aptima® Multitest Swab Specimen Collection Kit or Aptima® Unisex Swab Specimen Collection Kit

Collection: Vaginal swab: Collect vaginal fluid sample using the Aptima® swab by contacting the swab to the lower third of the vaginal wall and rotating the swab for 10 to 30 seconds to absorb fluid. Immediately place the swab into the transport tube and carefully break the swab shaft against the side of the tube. Tightly screw on the cap.

Causes for Rejection: Specimen with incorrect patient identification; unlabeled specimen; inappropriate specimen transport conditions, including specimens received frozen; specimens received after prolonged delay (usually >72 hours); specimen leaked in transit; specimen in expired transport or incorrect transport device; specimens with inappropriate source for test requested; specimen with fixative or additives; Aptima® urine transport; Aptima® swab transport >30 days from collection; Aptima® swab specimen without a swab; cleaning swab (white-shaft swab) in Aptima® swab transport; any non-Aptima® swab submitted in Aptima® transport device; transport device with multiple swabs; bloody or grossly mucoid specimens; bacterial swabs; specimen in ProbeTecTM UPT transport; ProbeTecTM Q-swabs

Storage Instructions: Maintain specimen at room temperature or refrigerate (2°C to 30°C).

Stability Requirements: Room temperature 30 days, Refrigerated 30 days, Frozen Not acceptable.

GENERAL INFORMATION

Expected TAT

3-4 days

Notes

RESULT CODE	RESULT NAME	LOINC
5195307	Atopobium vag	69565-0
5195308	BVAB 2	69566-8
5195309	Megasphaera 1	69567-6

CPT Code(s) 87798x3

Service Provided By

