

**STI Profile, CT/NG**


Order Name: **STI Profile, CT/NG**  
Test Number: 5195301  
Revision Date: 12/06/2023

TEST NAME	METHODOLOGY	LOINC CODE
STI Profile, CT/NG	<u>Assay Dependant</u>	

**SPECIMEN REQUIREMENTS**

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	See Instructions	See Instructions	See Instructions	Room Temperature
Instructions	<p><b>Specimen:</b> Serum and Swab or Urine <b>Volume:</b> Swab or 2 mL urine, 8 mL serum <b>Container:</b> Aptima® (Gen-Probe®) Endocx/ureth swab or urine, gel-barrier tube</p> <p><b>Collection:</b> Serum specimen: Usual blood collection technique. Ct/GC: Option 1: Gen-Probe® Aptima® endocervical or vaginal swab. Endocervical swab: Remove excess mucus from the cervical os and surrounding mucosa using the cleaning swab (white-shaft swab in the package with red printing). Discard this swab. Insert the specimen collection swab (blue-shaft swab in the package with green printing) into the endocervical canal. Gently rotate the swab clockwise for 10 to 30 seconds in the endocervical canal to ensure adequate sampling. Withdraw the swab carefully; avoid contact with the vaginal mucosa. Remove the cap from the swab specimen transport tube and immediately place the specimen collection swab into the transport tube. Carefully break the swab shaft at the scoreline; use care to avoid splashing of the contents. Recap the swab specimen transport tube tightly. Vaginal Swab: Care Provider Specimen: Collect vaginal fluid sample using the Gen-Probe® Aptima® vaginal swab kit by contacting the swab to the lower third of the vaginal wall and rotating the swab for 10 to 30 seconds to absorb fluid. Immediately place the swab into the transport tube and carefully break the swab shaft against the side of the tube. Tightly screw on the cap. Patient Self-Collection: Partially open the package of the Gen-Probe® Aptima® vaginal swab kit. Do not touch the soft tip or lay the swab down. If the soft tip is touched, the swab is laid down, or the swab is dropped, use a new Aptima® swab specimen collection kit. Remove the swab. Carefully insert the swab into the vagina about two inches past the introitus and gently rotate the swab for 10 to 30 seconds, making sure the swab touches the walls of the vagina so that moisture is absorbed by the swab. Withdraw the swab without touching the skin. Immediately place the swab into the transport tube, and carefully break the swab shaft against the side of the tube. Tightly screw on the cap. Ct/GC: Option 2: Urine Specimen: The patient should not have urinated for at least one hour prior to specimen collection. Direct patient to provide a first-catch urine (approximately 20 to 30 mL of the initial urine stream) into a urine collection cup free of any preservatives. Collection of larger volumes of urine may result in specimen dilution that may reduce test sensitivity; lesser volumes may not adequately rinse organisms into the specimen. Female patients should not cleanse the labial area prior to providing the specimen. Add urine to the Aptima® Combo 2 urine collection device. The final volume must be between the two black lines on the device (about 2 mL).</p>			

**GENERAL INFORMATION**

Expected TAT	3-5 Days
Notes	<b>Test Includes:</b> RPR with reflex to quantitative RPR and confirmatory Treponema pallidum antibodies; HIV p24 Antigen/Antibody with reflex to confirmation; HBV screening and diagnosis (HBsAg, anti-HBs, HBc-total with reflex to IgM); HCV antibody with reflex to quantitative real-time PCR; Chlamydia trachomatis NAA and Neisseria gonorrhoeae NAA. If reflex testing is performed, additional charges/CPT code(s) may apply.
CPT Code(s)	86592; 86704; 86706; 86803; 87340; 87389; 87491; 87591
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.