


Inhibitor Screen

Order Name: **INHIB SCR**
Test Number: 1516550
Revision Date: 09/09/2025

| TEST NAME | METHODOLOGY | LOINC CODE |
|---------------------------------|----------------|------------|
| PT 1:1 Mixture (Room Temp) | Clot Detection | |
| PTT 1:1 Mixture (Room Temp) | Clot Detection | |
| PTT-LA 1:1 Mixture (Room Temp) | Clot Detection | 67097-6 |
| PT 1:1 Mixture (Incubated) | Clot Detection | |
| PTT 1:1 Mixture (Incubated) | Clot Detection | |
| PTT-LA 1:1 Mixture (Incubated) | Clot Detection | |
| Inhibitor Screen Interpretation | | |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|--|--------------------|--------------------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 22 mL (16 mL) | Whole Blood | Sodium Citrate 3.2% (Blue Top) | Room Temperature |
| Alternate 1 | 12 mL (9 mL) | Double Spun Plasma | Sterile, Capped Plastic Tube | Frozen |
| Instructions | <p>Please list the patient's anticoagulant on the "Coagulopathy Questionnaire Form" and submit with specimen or fax to 918-744-2897.</p> <p>Please indicate anticoagulant therapy.</p> <p>Collect 6-8 Sodium Citrate 3.2% (Blue Top) tubes.</p> <p>Each 2.7mL Sodium Citrate 3.2% (Blue Top) tube must be filled to the proper level, no hemolysis. Improperly filled tubes can give erroneous results.</p> <p>Whole blood must be transported to lab immediately.</p> <p>If testing cannot be started within 4 hours of collection the specimen must be double spun then 1.5mL plasma aliquot from each tube into individual plastic aliquot tubes and freeze.</p> <p>Do not pool aliquots together!</p> <p>Coagulopathy Questionnaire Form</p> <p>Double Spin Procedure</p> | | | |

| GENERAL INFORMATION | |
|---------------------|--|
| Testing Schedule | Mon - Wed, Fri - Day Shift |
| Expected TAT | 2-4 Days |
| Clinical Use | Specific factor inhibitors are immunoglobulins with specificity for a single coagulation protein. The most common specific inhibitors are antibodies produced in relation to factor VIII. Nonspecific inhibitors, such as lupus anticoagulants, are also detected. |
| CPT Code(s) | If PT Mix is performed 85611X2 If PTT Mix is performed 85732X2 If PTT-LA Mix is performed 85732X2 Testing includes a pathology interpretation. |
| Service Provided By |  labcorp Oklahoma, Inc. |



**Ascension
St. John**

St. John Health System
Lab Catalog