


Creatinine

Order Name: **CREATININE**
Test Number: 2025050
Revision Date: 09/17/2025

TEST NAME	METHODOLOGY	LOINC CODE
Creatinine	<u>Kinetic Alkaline Picrate (Jaffe)</u>	2160-0

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated
Alternate 1	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated
Instructions	Specimen Stability: Ambient 7 days, Refrigerated 7 days, Frozen 3 months			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful for sensitive and quantitative measurement of renal function. A Glomerular Filtration Rate (GFR) is also provided. Reference range for GFR African American and non-African American Normal GFR: <60 mL/Min. Moderately decreased GFR: 30-59 mL/Min. Severely decreased GFR: 15-29 mL/Min. Kidney failure (or Dialysis): < 15 mL/Min.
Notes	Serum creatinine values will be used to calculate GFR (Glomerular Filtration Rate) with the 2021 CKD-EPI creatinine GFR calculation. The following test will be added on and reported at no additional charge. Order code: GFR Calculation [6907615] Result code: GFR CALC [6907615]
CPT Code(s)	82565
Service Provided By	 labcorp Oklahoma, Inc.