

## Von Willebrand Factor Antigen

Order Name: **VON WIL AG**  
Test Number: 5502200  
Revision Date: 10/08/2024

TEST NAME	METHODOLOGY	LOINC CODE
Von Willebrand Factor Antigen	<a href="#">Immunoturbidimetry</a>	27816-8

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (1.5mL)	Plasma	Sodium Citrate 3.2% (Blue Top)	Frozen
Instructions	Tubes must be filled to the proper level, no hemolysis. Improperly filled tubes can give erroneous results. Whole blood must be transported to lab immediately. <b>If testing cannot be started within 4 hours of collection</b> the specimen must be double spun to create One to Two 1.5mL plasma aliquot from each tube into individual plastic aliquot tubes and freeze. <b>Do not pool aliquots together!</b> Do not thaw. Hemolyzed specimens are not acceptable. <b>Overnight fasting is preferred.</b>			

GENERAL INFORMATION	
Testing Schedule	Thursday - Day Shift
Expected TAT	1-7 Days
CPT Code(s)	85246
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.