

**Tocainide Level**

Order Name: **TOCAINIDE**  
Test Number: 3611150  
Revision Date: 12/12/2022

| TEST NAME       | METHODOLOGY   | LOINC CODE |
|-----------------|---|------------|
| Tocainide Level | <u>Liquid Chromatography/Tandem Mass Spectrometry</u> |            |

| SPECIMEN REQUIREMENTS |   |               |   |                       |
|-----------------------|---|---------------|---|-----------------------|
| Specimen              | Specimen Volume (mL)  | Specimen Type | Specimen Container                        | Transport Environment |
| Preferred             | 1 mL (0.5 mL)   | Serum         | Clot Activator (Red Top, No-Gel)          | Room Temperature      |
| Alternate 1           | 1 mL (0.5 mL)   | Plasma        | Lithium Heparin (Dark Green Top / No-Gel) | Room Temperature      |
| Alternate 2           | 1 mL (0.5 mL)   | Plasma        | EDTA (Lavender Top)                       | Room Temperature      |
| Instructions          | <p><b>Specimen Type:</b> Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended.</p> <p><b>Specimen Storage:</b> SUBMISSION/TRANSPORT (&lt;3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.</p> <p><b>Specimen Collection:</b> Serum or plasma should be separated from cells within two hours of venipuncture. Submit serum or plasma in a plastic transport tube.</p> <p><b>Special Instructions:</b> Testing referred to MedTox Laboratories Inc tc 121</p> |               |   |                       |

| GENERAL INFORMATION |  |
|---------------------|--|
| Expected TAT        | 5 - 7 days   |
| Notes               | Labcorp Test Code: 811372  |
| CPT Code(s)         | 80299  |
| Service Provided By |  <b>labcorp</b><br>Oklahoma, Inc. |