

Pyruvic Acid, CSF

Order Name: **PYRU CSF**
Test Number: **6905131**
Revision Date: **12/12/2022**

| TEST NAME | METHODOLOGY | LOINC CODE |
|-------------------|-------------------------------|------------|
| Pyruvic Acid, CSF | <u>Quantitative Enzymatic</u> | 14122-6 |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|--|---------------------------|----------------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1 mL (0.5) | CSF (Cerebrospinal Fluid) | Sterile Screwtop Container | Frozen |
| Instructions | <p>Specimen Type: Sterile screw capped vial</p> <p>Specimen Storage: Freeze as soon as possible after collection</p> <p>Specimen Collection: Collect entire sample into a single sterile tube. Immediately after collection, add a 2:1 ratio (acid:sample) of COLD 8% perchloric acid. Vortex for 30 seconds, and place on ice for 5 minutes. Store frozen at -20 degrees Celsius and ship frozen.</p> <p>Special Instructions: Testing referred to Medical Neurogenetic LLC MNEGA#MET11</p> <p>Specimen Stability: Ambient: unstable (stability provided by manufacturer or li, Refrigerated : 24 hours (stability provided by manufacturer or li, Frozen: -20 C = 72 hours; -80 C = Indefinite (stability pr</p> | | | |

| GENERAL INFORMATION | |
|---------------------|--|
| Expected TAT | 10 - 14 days |
| Notes | Labcorp Test Code: 620045 |
| CPT Code(s) | 84210 |
| Service Provided By |  labcorp Oklahoma, Inc. |