

**Triiodothyronine, Free by Equilibrium Dialysis/HPLC-Tandem Mass Spectrometry**

Order Name: **T3 FREE SP**  
Test Number: 6905173  
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Triiodothyronine, Free by Equilibrium Dialysis/HPLC-Tandem Mass Spectrometry	<u>Quant Equilibrium Dialysis/HPLC-Tandem Mass Spectrometry</u>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.3 mL)	Serum	Clot Activator (SST or Red No-Gel)	Frozen
Alternate 1	1 mL (0.3 mL)	Plasma	EDTA (Lavender Top)	Frozen
Alternate 2	1 mL (0.3 mL)	Plasma	Sodium Heparin (Green Top, No-Gel)	Frozen
Instructions	<p><b>Notes:</b> 0.3 mL (Note: This volume Does NOT allow for repeat testing).</p> <p><b>Specimen Type:</b> Red-top tube, gel-barrier tube, lavender-top (EDTA) tube or green-top (sodium heparin) tube</p> <p><b>Specimen Storage:</b> <b>FREEZE</b> (preferred) or refrigerate Stability: Room temperature up to 5 days Refrigerated up to 5 days Frozen: up to 70 days</p> <p><b>Specimen Collection:</b> <b>Serum/plasma must be separated from cells within 45 minutes of venipuncture, transfer into a plastic transfer tube and Freeze.</b></p> <p><b>Special Instructions:</b> Testing is performed at Esoterix Endocrinology Laboratory: ES#500336</p> <p><b>Specimen Stability:</b> Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available</p>			

GENERAL INFORMATION	
Expected TAT	6 - 10 days
Notes	Labcorp Test Code: 503600
CPT Code(s)	84481
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.