


**Vibrio Culture**

Order Name: **C VIBRIO**  
Test Number: 6001076  
Revision Date: 07/10/2017

TEST NAME	METHODOLOGY	LOINC CODE
Vibrio Culture	Culture	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3mL	Stool, Random	ETM or Cary Blair container	Ambient / Refrigerated
Instructions	Stool specimens should be placed in Cary-Blair based transport media immediately after collection (ETM or equivalent). Add enough stool to fill container to indicator line. Specimen Stability: Room Temperature: 72 hours, Refrigerated: 72 hours, Frozen: Unacceptable Unacceptable Specimens: Stool specimens must not be submitted in formalin or PVA preservative and should not contain barium or oily cathartics. Specimens outside of appropriate transport media for 4 or more hours are unacceptable. Stools received 72 hours after collection or specimens that have been frozen are unacceptable. Multiple specimens collected within 24 hours will not be tested.			

GENERAL INFORMATION	
Testing Schedule	Sun - Sat
Expected TAT	2-3 Days
CPT Code(s)	87046; Additional CPT codes if identification is required.
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.