


**Complement, Total (CH50)**

Order Name: **CH50**  
Test Number: **5569251**  
Revision Date: **04/14/2023**

TEST NAME	METHODOLOGY	LOINC CODE
Complement, Total (CH50)	<u>Quantitative IVD Assay</u>	4532-8

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.3 mL)	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Instructions	<p>Allow 1 hour to clot at room temperature, Separate the serum from cells ASAP or no longer than 2 hours after collection.</p> <p>Separate the serum from the cells and transfer 1mL(0.3mL) of serum into a standard transport tube and freeze within 2 hours of collection.<b>Do Not Share this specimen with another test.</b></p> <p><b>It is critical to freeze the complement specimen immediately after the transfer of specimen to a transport tube.</b></p> <p>Reasons for Rejection: Use of serum separator tube, Clotting at 2-8°C, Exposure to repeated freeze/thaw cycles, Samples containing high levels of lipid, hemoglobin or bilirubin cause interference and should be avoided.</p> <p>Stability: Room Temperature 2 hours, Refrigerated n/a, Frozen 2 Weeks</p>			

GENERAL INFORMATION	
Testing Schedule	Tuesday, Friday
Expected TAT	5-7 days, batched
Clinical Use	Low levels of total complement may occur during infection, exacerbation of SLE, exacerbation of hereditary angioedema and glomerulonephritis. Undetectable levels suggest possibility of a complement deficiency.
Notes	This assay is an FDA approved quantitative IVD assay that will improve testing reliability, workflow and turn around time (TAT). The methodology utilizes liposomes, encapsulating glucose-6-phosphate dehydrogenase (G6PDH) to mimic an invading microorganism. The patient's serum sample provides the complement for the lysis of the liposome which releases G6PDH and reacts with glucose-6-phosphate and NAD. The change in absorbance is measured and is proportional to the complement activity in the sample. The results are compared to a calibration curve which gives a value for the patient's sample. Complement activity has been correlated with the active stage of systemic lupus erythematosus, rheumatoid arthritis, some forms of nephritis and inherited deficiencies of the classical complement system.
CPT Code(s)	86162
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.