


**HPV High Risk with 16/18 Genotype - SurePath**

Order Name: **HPV HI/16/18 SP**  
Test Number: **5522580**  
Revision Date: **01/01/2025**

TEST NAME	METHODOLOGY	LOINC CODE
HPV Genotype 16 PCR	Qualitative PCR	61372-9
HPV Genotype 18 PCR	Qualitative PCR	61373-7
HPV Other High Risk	Qualitative PCR	70061-7
HPV Source	Prompt	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (1 mL)	PAP specimen	ThinPrep PreservCyt solution	Room Temperature
Alternate 1	2 mL (1 mL)	PAP specimen	SurePath Liquid Pap Container (Pap Prep)	Room Temperature
Instructions	<p><b>Notes:</b> 1 mL (Note: This volume Does NOT allow for repeat testing.)</p> <p><b>Specimen Type:</b> ThinPrep(R) vial or SurePath(TM) vial</p> <p><b>Specimen Storage:</b> Maintain liquid-based cytology specimens at room temperature. date of collection prior to performing the cobas(R) HPV test. ThinPrep(R) specimens should not be frozen.</p> <p><b>Specimen Collection:</b> BRUSH/SPATULA TECHNIQUE: Insert the brush into the endocervical canal until only the bottommost fibers are exposed. Slowly rotate the brush 1/4 to 1/2 turn in one direction. Do NOT over-rotate the brush. Then, rotate the brush in the PreservCyt(R) solution 10 times while pushing against the wall of the ThinPrep(R) vial. Swirl the brush vigorously to release additional material. Discard the brush. Obtain an adequate sample from the ectocervix using a plastic spatula. Swirl vigorously in the ThinPrep(R) vial 10 times and discard the spatula. Tighten the cap on the ThinPrep(R) container so that the torque line on the cap passes the torque line on the vial. SUREPATH(TM) VIAL: When using the SurePath(TM) vial, the cervical broom must be used for specimen collection. Insert the broom into the cervical os and rotate five times. Place the broom head into the CytoRich(TM) preservative fluid in the SurePath(TM) collection vial. Tightly cap the vial.</p> <p><b>Specimen Stability:</b> Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available</p>			

GENERAL INFORMATION	
Expected TAT	2 - 4 days
Clinical Use	<p>This test specifically identifies types HPV16 and HPV18 while concurrently detecting the rest of the high-risk types: 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68 without further specific differentiation.</p> <p>Limitations: Detection of high-risk HPV is dependent on the number of copies present in the specimen and may be affected by specimen collection methods, patient factors, stage of infection, and the presence of interfering substances.</p>
CPT Code(s)	87626
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.