

Ct/GC/Tv NAA+M genitalium Urine

Order Name: **Ct/GC/Tv M.gen Ur**
 Test Number: 5195318
 Revision Date: 05/01/2026

TEST NAME	METHODOLOGY	LOINC CODE
Ct/GC/Tv NAA+M genitalium Urine	<u>Nucleic acid amplification (NAA)</u>	

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	30 mL (5 mL)	Urine (See Inst.)	APTIMA Urine Transport Tube	Room Temperature

Instructions

Specimen: First-void urine in Aptima® urine transport
Volume: 5 mL or 30 mL urine collection
Container: Aptima® urine transport
Collection: Urine specimen: The patient should not have urinated for at least one hour prior to specimen collection. Direct patient to provide a first-catch urine (approximately 20 to 30 mL of the initial urine stream) into a urine collection cup free of any preservatives. Collection of larger volumes of urine may result in specimen dilution that may reduce test sensitivity; lesser volumes may not adequately rinse organisms into the specimen. Female patients should not cleanse the labial area prior to providing the specimen. Add urine to the Aptima® Combo 2 urine collection device. The final volume must be between the two black lines on the device (about 2 mL).
Storage Instructions: Maintain specimen at room temperature or refrigerate (2°C to 30°C).
Stability Room temperature 30 days (stability provided by manufacturer or literature reference)

GENERAL INFORMATION

Expected TAT 3-4 days

Notes	RESULT CODE	RESULT NAME	LOINC CODE
	6907667	Mycoplasma genitalium Urine	88226-6
	6910159	Trichomonas Urine	70167-2
	6910151	Chlamydia Urine	42931-6
	6910155	Gonorrhea Urine	60256-5

Service Provided By

