

Eosinophil Count, Absolute

Order Name: **EOS CT ABS**Test Number: 0100050
Revision Date: 10/31/2019

TEST NAME			METHODOLOGY	LOINC CODE	
Eosinophil Count, Absolute			Flow cytometry	711-2	
SPECIMEN REQU	UREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	5 mL (1 mL)	Whole Blood	EDTA (Lavender Top)	Refrigerated	
Alternate 1	1 mL (0.5 mL)	Whole Blood	EDTA (Lavender) Microtainer/Bullet	Refrigerated	
Instructions	For best results				
	Room temperature specimens should be tested within 24hrs, otherwise send Refrigerated. Refrigerated specimens can be tested up to 48hrs.				
Specimens received greater than 48hrs old will be canceled.					

GENERAL INFORMATION		
Testing Schedule	Daily	
Expected TAT	1 Day	
Clinical Use	Eosinophilia is found in parasitic diseases, allergic diseases, skin disorders, and certain malignancies.	
CPT Code(s)	85004	
Service Provided By	labcorp Oklahoma, Inc.	