


Tyrosine Urine, Random

Order Name: **TYROSINE U**  
Test Number: 1002000  
Revision Date: 08/26/2014

TEST NAME	METHODOLOGY	LOINC CODE
Tyrosine Urine, Random	Visual *****	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (1 mL)	Urine, Random	Sterile Urine container	Refrigerated
Instructions	Early morning specimens preferred. Room temperature is acceptable is tested immediately, otherwise Transport or Store the specimen Refrigerated.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1 Day
CPT Code(s)	84510
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.