

Anti-Thrombin 3 (ATIII), Functional Activity

Order Name: **THROM3 FUN**Test Number: 1501825
Revision Date: 05/14/2025

TEST NAME			DDOLOGY	LOINC CODE	
Anti-Thrombin 3 (ATIII), Functional Activity		Chrom	ogenic	27811-9	
SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	2.7 mL	Whole Blood	Sodium Citrate 3.2% (Blue Top)	Room Temperature	
Alternate 1	1.5 mL	Double Spun Plasma	Sodium Citrate 3.2% (Blue Top)	Frozen	
Instructions	Please indicate anticoagulant therapy. Each 2.7mL Sodium Citrate 3.2% (Blue Top) tube must be filled to the proper level, no hemolysis. Improperly filled tubes can give erroneous results. Whole blood must be transported to lab immediately. If testing cannot be started within 4 hours of collection the specimen must be double spun then 1.5 mL plasma aliquot from each tube into individual plastic aliquot tubes and freeze. Do not pool aliquots together! Do not thaw.				

GENERAL INFORMATION		
Testing Schedule	Tues, Thurs	
Expected TAT	1-7 Days	
Clinical Use	Antithrombin III is used to assess thrombotic risk.	
CPT Code(s)	85300	
Service Provided By	Oklahoma, Inc.	