


Complement C4, Serum

Order Name: **C4**
Test Number: 5000350
Revision Date: 02/26/2018

TEST NAME	METHODOLOGY	LOINC CODE
Complement C4, Serum	<u>Turbidometric</u>	4498-2

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	2mL (1mL)	Serum	Clot Activator SST	Refrigerated or Frozen
Alternate 1	2mL (1mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated or Frozen
Instructions	Separate serum from cells within one hour of collection. Stability: Room Temperature N/A, Refrigerated 7 days, Frozen 8 days. If ordering C4 with CH50, Please Collect a Red No-Gel Clot tube, freeze two separate individual aliquots for C4 and CH50. (avoid repeated freeze/thaw cycles).			

GENERAL INFORMATION	
Testing Schedule	Mon-Sat
Expected TAT	3 Days
Clinical Use	Low levels may occur during infections, exacerbation of SLE and glomerulonephritis. Undetectable level suggests C4 deficiency.
CPT Code(s)	86160
Service Provided By	 labcorp Oklahoma, Inc.