


Haptoglobin

Order Name: **HAPTOGLOBN**
Test Number: **5000700**
Revision Date: **05/20/2021**

TEST NAME	METHODOLOGY	LOINC CODE
Haptoglobin	<u>Immunoturbidimetry</u>	4542-7

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated
Alternate 1	1 mL (0.5)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated
Instructions	Stability: Ambient 4 hours, Refrigerated 7 days, Frozen 2 weeks			

GENERAL INFORMATION	
Testing Schedule	Mon - Fri
Expected TAT	1-3 days
Clinical Use	A low serum haptoglobin concentration is not specific for hemolysis. Evaluation and diagnosis of hemolytic anemia.
CPT Code(s)	83010
Service Provided By	 labcorp Oklahoma, Inc.