

**Anti-Nuclear Antibody Screen (ANA) with Reflex to Titer**


Order Name: **ANA SCR**  
Test Number: 5500050  
Revision Date: 11/12/2020

| TEST NAME   | METHODOLOGY                          | LOINC CODE |
|---|--------------------------------------|------------|
| Anti-Nuclear Antibody Screen (ANA) with Reflex to Titer | <u>Indirect Fluorescent Antibody</u> | 8061-4     |

| SPECIMEN REQUIREMENTS |  |               |                                    |                       |
|-----------------------|--|---------------|------------------------------------|-----------------------|
| Specimen              | Specimen Volume (min)  | Specimen Type | Specimen Container                 | Transport Environment |
| Preferred             | 1 mL   | Serum         | Clot Activator (SST or Red No-Gel) | Refrigerated          |
| Instructions          | Allow specimen to clot completely at room temperature. <b>Separate serum from cells ASAP or within 2 hours of collection.</b><br>Stability After separation from cells: Ambient 8 hours, Refrigerated 7 days, Frozen 3 months (avoid repeated freeze/thaw cycles). |               |                                    |                       |

| GENERAL INFORMATION |  |
|---------------------|--|
| Testing Schedule    | Mon - Fri  |
| Expected TAT        | 3 Days   |
| Clinical Use        | To aid in the diagnose autoimmune diseases.  |
| Notes               | If ANA screen is positive >1:80; titer is automatically performed. ANA screening is done by IFA using HEp-2 Cells.<br>See <a href="https://www.anapatterns.org/index.php">https://www.anapatterns.org/index.php</a> for current clinical associations and identified patterns. |

| REFLEX ORDER CODE | REFLEX TEST NAME           | RESULT CODE | RESULT LONG NAME               |
|-------------------|----------------------------|-------------|--------------------------------|
| 5512450           | Antinuclear Antibody Titer | 5500100     | Antinuclear Antibody Titer     |
|                   |                            | 5500150     | Antinuclear Antibody Pattern   |
|                   |                            | 5577550     | Antinuclear Antibody Titer 2   |
|                   |                            | 5577650     | Antinuclear Antibody Pattern 2 |
|                   |                            | 5577600     | Antinuclear Antibody Titer 3   |
|                   |                            | 5577700     | Antinuclear Antibody Pattern 3 |

|                     |  |
|---------------------|--|
| CPT Code(s)         | 86038 (86039 will bill in addition if reflexed to titer)   |
| Service Provided By |  <b>labcorp</b><br>Oklahoma, Inc. |