

**Smooth Muscle Antibody Titer**

Order Name: **SM MUSC DE**  
Test Number: **5572400**  
Revision Date: **10/23/2017**

TEST NAME	METHODOLOGY	LOINC CODE
Smooth Muscle Antibody Titer	<u>Indirect Fluorescent Antibody</u>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1.0 mL	Serum	Clot Activator SST	Refrigerated
Instructions	Allow specimen to clot completely at room temperature. <b>Separate serum or plasma from cells ASAP or within 2 hours of collection.</b> Stability After separation from cells: Ambient 4hours, Refrigerated 7 days, Frozen 1 month (avoid repeated freeze/thaw cycles).			

GENERAL INFORMATION	
Testing Schedule	Mon - Fri
Expected TAT	3 Days
Clinical Use	Smooth muscle antibody is in high titers ( $\geq$ 1:160) in approximately 97% of patients with autoimmune chronic active hepatitis.
CPT Code(s)	86256
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.