


**Heparin Induced Platelet Antibody (HIT)**Order Name: **HEPARIN AB**Test Number: **5566580**Revision Date: **08/19/2020**

TEST NAME	METHODOLOGY	LOINC CODE
Heparin Induced Platelet Antibody	<u>Quantitative Enzyme-Linked Immunosorbent Assay</u>	57761-9
Optical Density		
% inhibition	<u>Calculation</u>	

**SPECIMEN REQUIREMENTS**

Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL	Serum	Clot Activator SST	Frozen
Instructions	Centrifuge, separate and freeze serum <b>Stability:</b> 2-8 C for up to 48 hours, or frozen			

**GENERAL INFORMATION**

Testing Schedule	Mon - Fri
Expected TAT	3 Days
Clinical Use	The heparin platelet antibody test is designed to detect antibodies against platelet factor 4 (PF4) that are created when PF4 is complexed with heparin for heparin induced thrombocytopenia (HIT).
CPT Code(s)	86022
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.