

Chlamydia Psittaci Antibody

Order Name: CHLAM PSIT
Test Number: 5571525
Revision Date: 06/09/2003

| TEST NAME | METHODOLOGY | LOINC CODE |
|-----------------------------------|-------------------------------|------------|
| Chlamydia psittaci IgG | Indirect Fluorescent Antibody | 43868-9 |
| Chlamydia psittaci IgM | Indirect Fluorescent Antibody | 44977-7 |
| Chlamydia Psittaci Interpretation | | 44081-8 |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|-----------------------|---------------|--------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1 mL | Serum | Clot Activator SST | Refrigerated |

| GENERAL INFORMATION | |
|---------------------|---|
| Testing Schedule | Mon - Fri |
| Expected TAT | 3 Days |
| Clinical Use | C. psittaci causes a respiratory illness which is usually found in pet bird owners or poultry industry owners. Interpret the serologic results in light of clinical history and presentation. |
| CPT Code(s) | 86631, 86632 |
| Service Provided By | Oklahoma, Inc. |