


Anti-Nuclear Antibody (ANA) Serous Fluid

Order Name: **ANA FL**
Test Number: **5590550**
Revision Date: **07/14/2005**

TEST NAME	METHODOLOGY	LOINC CODE
Anti-Nuclear Antibody (ANA) Serous Fluid	<u>Indirect Fluorescent Antibody</u>	14607-6

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL	Serous Fluid	Clot Activator (Red Top, No-Gel)	Room Temperature
Instructions	1 mL Serous Fluid			

GENERAL INFORMATION	
Testing Schedule	Mon - Fri
Expected TAT	3 Days
Notes	ANA screening is done by IFA using HEp-2 Cells.
CPT Code(s)	86038
Service Provided By	 labcorp Oklahoma, Inc.