


Gliadin Deamidated Antibody, IgA and IgG

Order Name: **GLIADIN A/G**
Test Number: 5558560
Revision Date: 01/01/2022

| TEST NAME | METHODOLOGY | LOINC CODE |
|-----------------------------|--------------------------------|------------|
| Gliadin Antibody IgA | <u>Chemiluminescence Assay</u> | 63453-5 |
| Gliadin Antibody IgG | <u>Chemiluminescence Assay</u> | 63459-2 |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|---|---------------|---------------------------|-----------------------|
| Specimen | Specimen Volume (mL) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 2mL (1mL) | Serum | Clot Activator SST | Frozen |
| Instructions | Allow specimen to clot completely at room temperature. Separate serum from cells ASAP or within 2 hours of collection. Stability: Room Temperature 8hrs, Refrigerated: 48hrs, Frozen: 30 days. Transport Frozen | | | |

| GENERAL INFORMATION | |
|----------------------------|--|
| Testing Schedule | Mon - Fri |
| Expected TAT | 1-3 Days |
| Clinical Use | Assist in diagnosis of celiac disease; however, interpret results with caution due to the propensity of assay false positives. Useful to monitor diet compliance in celiac patients. |
| CPT Code(s) | 86258x2 |
| Service Provided By |  labcorp Oklahoma, Inc. |