


Ketone Bodies Qualitative

Order Name: **KETONES**
Test Number: 2000075
Revision Date: 03/23/2021

TEST NAME	METHODOLOGY	LOINC CODE
Ketone Bodies Qualitative	<u>Dry Chemistry</u>	2513-0

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated
Alternate 1	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated
Instructions	Keep tightly stoppered - Do not remove stopper of collection tube. Keep refrigerated. This test is only available on inpatients at St. John Sapulpa, St. John Owasso, St. John Broken Arrow and Jane Phillips Medical Center. All other locations please use Beta Hydroxybutyrate .			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful for the diagnosis and monitoring of therapy for diabetic ketoacidosis.
CPT Code(s)	82009
Service Provided By	 labcorp Oklahoma, Inc.