


Bilirubin, Total And Direct

Order Name: **BILI T/D**  
Test Number: 2001000  
Revision Date: 04/26/2021

| TEST NAME               | METHODOLOGY     | LOINC CODE |
|-------------------------|-----------------|------------|
| <b>Bilirubin Direct</b> | Diazo           | 1968-7     |
| <b>Bilirubin Total</b>  | Jendrassik-Grof | 1975-2     |

| SPECIMEN REQUIREMENTS |   |               |                                       |                       |
|-----------------------|---|---------------|---------------------------------------|-----------------------|
| Specimen              | Specimen Volume (min)   | Specimen Type | Specimen Container                    | Transport Environment |
| Preferred             | 1 mL (0.5)  | Plasma        | Lithium Heparin PST (Light Green Top) | Refrigerated          |
| Alternate 1           | 1 mL (0.5)  | Serum         | Clot Activator SST                    | Refrigerated          |
| Instructions          | Exposure to direct sunlight can decrease bilirubin in samples up to 50% within one hour.<br>Light protected preferred; specimens not protected from light but refrigerated/frozen will be accepted.<br>Stability: Ambient 1 day (if protected from light), Refrigerated 7 days, Frozen 3 months |               |                                       |                       |

| GENERAL INFORMATION |  |
|---------------------|--|
| Testing Schedule    | Daily  |
| Expected TAT        | 1-2 days   |
| Clinical Use        | Useful in evaluating hepatocellular diseases, hepatitis, cirrhosis, and jaundice.                                    |
| CPT Code(s)         | 82248 - Bilirubin direct 82247 - Bilirubin total   |
| Service Provided By |  <b>labcorp</b><br>Oklahoma, Inc. |