

Order Name: BUN POST
Test Number: 2009750
Revision Date: 06/04/2021

TEST NAME			METHODOLOGY	LOINC CODE
Urea Nitrogen, Blood Post Dialysis		Urease/GLDH	11064-3	
SPECIMEN REQU	IREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated
Alternate 1	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated
Instructions	Stability: Ambient 7 days, Refu	igerated 7 days, Frozen 1 ye	ear	

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Post dialysis. Useful in the monitoring of dialysis.
CPT Code(s)	84520
Service Provided By	labcorp Oklahoma, Inc.