


**Cyclosporine Level Peak**

Order Name: **CYCLO PEAK**  
Test Number: 2015325  
Revision Date: 05/03/2021

TEST NAME	METHODOLOGY	LOINC CODE
Cyclosporine Level Peak	<u>Chemiluminescent Micropartical Immunoassay</u>	53834-8

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (1.0)	Whole Blood	EDTA (Lavender Top)	Refrigerated
Instructions	Peak: draw 2 hours after medication. Stability: Ambient 24 hours. Refrigerated 7 days. Frozen > 7 days.			

GENERAL INFORMATION	
Testing Schedule	Mon - Fri
Expected TAT	1-2 days
Clinical Use	This test is used to monitor blood concentration of Cyclosporine A.
CPT Code(s)	80158
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.